


MAJESTROUGH[®]

Peripheral Interventional



DOC

Diverse

Lubricity

Seamless


Stability

let us guide you. . .

Hicicare Science Co., Ltd.

☎ 0086-20-84969556

Structural Components



Product Features

- 1 Meets various clinical needs
- 2 Hydrophilic coating
- 3 1:1 Torque Control
- 4 Best-in-Class Visibility
- 5 Smooth plasticity
- 6 DOC extension

Ordering Information

GUIDEWIRE SERIES	SERIES NAME	FUNCTIONAL	TIP LOAD(g)	DOC
0.014"	SPEEDTRACK 14	Working guidewire	2.8	•
	SWIFT 14	Working guidewire	3.1	•
	AERO 14	Working guidewire	2.7	•
	SPEEDTRACK 14 ES	Support guidewire	3.5	•
	SPEEDTRACK 14 CTO	CTO-dedicated	2.7	•
0.016"	SUREPASS 16	Working guidewire	5.5	•
	SPEEDTRACK 16	Working guidewire	3.2	•
	AERO 16	Working guidewire	2.8	•
0.018"	SPEEDTRACK 18	Working guidewire	4	•
	SUREPASS 18	Working guidewire	6	•
	SPEEDTRACK FLEX	High-penetration	14.9	•
	SPEEDTRACK 18 CTO	CTO-dedicated	30	•