

PRIMTHROUGH[®]

Leading Guide Wire



Lubricity

Trackability

Tip Durability

Steerability


Stability

let us guide you. . .

Hicicare Science Co., Ltd.

☎ 0086-20-84969556

Structural Components



Product Features

- 1 Extra Floppy and Atraumatic tip
- 2 Hydrophilic coating
- 3 Reliable Rail Support
- 4 Best-in-Class Visibility
- 5 1:1 Torque Control for Steerability, Pushability and Trackability

Ordering Information

CATALOG NO.	LENGTH (CM)	TIP SHAPE	DIAMETER (INCH/MM)	COMPONENTS
GWN185C06S	185	Straight	0.014/0.36	—
GWN185C06J	185	J-Shape	0.014/0.36	—
GWN300C06S	300	Straight	0.014/0.36	—
GWN300C06J	300	J-Shape	0.014/0.36	—
GWN185C10S	185	Straight	0.014/0.36	—
GWN185C10J	185	J-Shape	0.014/0.36	—
GWN300C10S	300	Straight	0.014/0.36	—
GWN300C10J	300	J-Shape	0.014/0.36	—
GWN185S06S	185	Straight	0.014/0.36	Shaping Tool
GWN185S06J	185	J-Shape	0.014/0.36	Shaping Tool
GWN300S06S	300	Straight	0.014/0.36	Shaping Tool
GWN300S06J	300	J-Shape	0.014/0.36	Shaping Tool
GWN185S10S	185	Straight	0.014/0.36	Shaping Tool
GWN185S10J	185	J-Shape	0.014/0.36	Shaping Tool
GWN300S10S	300	Straight	0.014/0.36	Shaping Tool
GWN300S10J	300	J-Shape	0.014/0.36	Shaping Tool